Mölnlycke®

Mepilex® EM

WITH SAFETAC® TECHNOLOGY



Manufacturer

 Mölnlycke Health Care AB  
Gamlestadvägen 3C, Box 13080  
SE-402 52 Göteborg, Sweden

[www.molnlycke.com](http://www.molnlycke.com)

## Mepilex® EM

Single use  
Usage unique

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Soft silicone foam dressing

## Product description

Mepilex EM is a highly conformable dressing that absorbs exudates, maintains a moist wound environment and minimises the risk for maceration. The dressing has a Safetac® wound contact layer that is a unique adhesive technology. It minimises pain to patients and trauma to wounds and the surrounding skin at dressing removal.

Mepilex EM consists of

- a soft silicone wound contact layer (Safetac)
- a thin, flexible absorbent pad of polyurethane foam
- an outer polyurethane film which is breathable but waterproof

Dressing material content:

Silicone, polyurethane

## Indications for use

Mepilex EM is designed for the management of a wide range of non to low exuding wounds such as leg and foot ulcers, pressure ulcers, partial thickness burns, radiation skin reactions and Epidermolysis Bullosa. Mepilex EM can also be used as a protection of compromised and/or fragile skin.

## Precautions

- Do not use on patients with known hypersensitivity to the dressing or its materials.
- Do not use together with oxidising agents such as hypochlorite solutions or hydrogen peroxide.
- If you see signs of clinical infection e.g. fever or if the wound or surrounding skin becomes red, warm or swollen, consult a health care professional for appropriate treatment.
- Do not reuse, if reused performance of the product may deteriorate and cross contamination may occur.
- Sterile. Do not use if sterile barrier is damaged or opened prior to use. Do not re-sterilise.

## Instructions for use

Mepilex EM can be used by laypersons under the supervision of health care professionals.

1. Cleanse the wound in accordance with clinical practice.
2. Dry the surrounding skin thoroughly.
3. Select an appropriate dressing size. The dressing should cover the dry surrounding skin by at least 1-2 cm for small sizes (sizes up to 12.5x12.5 cm) and 3-5 cm for large sizes. If required, the dressing may be cut to suit various wound shapes and locations.
4. Remove the first release film and apply the adherent side to the wound.
5. Remove the remaining release film and smooth down the dressing on the skin. Do not stretch the dressing.
6. When necessary, fixate the dressing with a bandage or other fixation.

The dressing change interval may be several days. Change the dressing before it is fully saturated, at signs of leakage, or as indicated by clinical practice.

Mepilex EM can be used under compression bandaging and in combination with gels.

Disposal should be handled according to local environmental procedures.

## Other information

The polyurethane foam used in the product may change colour to more yellow when it is exposed to light, air and/or heat. The colour change has no influence on product properties when used before expiry date.

If any serious incident has occurred in relation to the use of Mepilex EM, it should be reported to Mölnlycke Health Care.

Mepilex® and Safetac® are registered trademarks of Mölnlycke Health Care AB.

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**CANADA** Mölnlycke Health Care Inc.,  
2010 Winston Park Drive, Suite 100,  
Oakville, Ontario, L6H 5R7  
Tel: +1 905 829 1502

**CHINA** Mölnlycke Health Care  
(Shanghai) Co., Ltd.  
RM 629, No.8 Hua Jing Road,  
Wai Gao Qiao FTZ Shanghai, China  
Tel: +86 10 5128 8571

**CZECH REPUBLIC** Mölnlycke Health Care s.r.o.  
Hajkova 2747/22, 130 00 Prague 3  
Tel: +420 221 890 517 (**Reception**)  
Tel: +420 221 890 511 (**Customer Service**)

**DENMARK** Mölnlycke Health Care ApS,  
Gydevang 33, DK-3450 Allerød  
Tel: +45 48 16 82 68  
**Customer Service** Tel: +80 886 810

**ESTONIA** Mölnlycke Health Care OÜ,  
Vabaõhumuuseumitee 4-67, EE-13522 Tallin  
Tel: +372 671 1520

**FINLAND** Mölnlycke Health Care Oy,  
Pitäjänmäentie 14, 3 krs, FI-00380 Helsinki  
Tel: +358 201 622 300  
**Customer Service** Tel: 0800 113 442

**FRANCE** Mölnlycke Health Care SAS,  
13 Allée du Château Blanc  
59290 WASQUEHAL  
Tel: +33 320 122 555  
**Customer Service** Tel: 0800 910 292

**GERMANY** Mölnlycke Health Care GmbH  
Grafenberger Allee 297  
DE-40237 Düsseldorf  
Tel: +49 211 92 08 80  
**Customer Service** Tel: 0800 186 21 80

**HUNGARY** Mölnlycke Health Care Kft  
HU-1134 Budapest, Dévai u. 26-28. III. em.  
Tel: +36 1 477 3080

**INDIA** Mölnlycke Health Care India Pvt Ltd  
Universal Business Centre  
Second Floor, Universal Trade Tower  
Sector 49, Gurgaon Sohna Road  
Gurgaon 122 018, Haryana  
Tel: +91 124 4696 008

**ITALY** Mölnlycke Health Care s.r.l.,  
Via Marsala 40/C, IT-21013 Gallarate (VA)  
Tel: +39 0331 714 011

**JAPAN** Mölnlycke Health Care KK,  
4F Concieria Nishi-Shinjuku Tower's West  
6-20-7 Nishi-Shinjuku,  
Shinjuku-Ku, Tokyo 160-0023  
Tel: +81 3 6914 5004

**KOREA** Mölnlycke Health Care Korea Co., Ltd  
2nd floor, Hansol Richville.  
46 Wiryeseong-daero,  
Songpa-gu, Seoul, 05627  
South Korea  
Tel: +82 2 3789 1402

**LATVIA** Mölnlycke Health Care AB,  
Dzelzavas 120 M, Riga LV-1021  
Tel: +372 6711 522

**LITHUANIA** Mölnlycke Health Care AB,  
J. Galvydzio g. 3, LT 08236 Vilnius  
Tel: +370 5 2274 58 37

**NETHERLANDS** Mölnlycke Health Care B.V.  
Postbus 3196, 4800 DD Breda  
Tel: +31 76 521 9663

**NORWAY** Mölnlycke Health Care AS  
Postboks 6229 Etterstad, NO-0603 Oslo  
Tel: +47 22 70 63 70  
**Customer Service** Tel: 800 161 78

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Mölnlycke Health Care AG
Mölnlycke Health Care Ltd.
**USA** Mölnlycke Health Care US, LLC